


Productlist 2008

Pins and Wires




CUTTING EDGE PINS: SELF DRILLING, SELF TAPPING




| Order Nr. | L= Length [mm] | A= Thread [mm] | Ø [mm] | S [mm] | AMOUNT | PRICE |
|-------------|-------------------|-------------------|-----------|-----------|--------|-------|
| F33- 363450 | 60 | 15 | 3,0 | 2,0 | 1 | |
| F33- 363451 | 60 | 20 | 3,0 | 2,0 | 1 | |
| F33- 363452 | 60 | 25 | 3,0 | 2,0 | 1 | |
| F33- 363454 | 80 | 20 | 3,0 | 2,0 | 1 | |
| F33- 363455 | 80 | 25 | 3,0 | 2,0 | 1 | |
| F33- 363456 | 80 | 30 | 3,0 | 2,0 | 1 | |
| F66- 364637 | 125 | 30 | 4,0 | 3,0 | 1 | |
| F66- 364638 | 125 | 35 | 4,0 | 3,0 | 1 | |
| F66- 364639 | 125 | 40 | 4,0 | 3,0 | 1 | |
| F66- 364413 | 150 | 30 | 4,0 | 3,0 | 1 | |
| F66- 364414 | 150 | 35 | 4,0 | 3,0 | 1 | |
| F66- 364415 | 150 | 40 | 4,0 | 3,0 | 1 | |
| F66- 365639 | 125 | 40 | 5,0 | 4,0 | 1 | |
| F66- 365640 | 125 | 45 | 5,0 | 4,0 | 1 | |
| F66- 365641 | 125 | 50 | 5,0 | 4,0 | 1 | |
| F66- 365415 | 150 | 40 | 5,0 | 4,0 | 1 | |
| F66- 365416 | 150 | 45 | 5,0 | 4,0 | 1 | |
| F66- 365417 | 150 | 50 | 5,0 | 4,0 | 1 | |
| F66- 365645 | 175 | 40 | 5,0 | 4,0 | 1 | |
| F66- 365646 | 175 | 45 | 5,0 | 4,0 | 1 | |
| F66- 365647 | 175 | 50 | 5,0 | 4,0 | 1 | |
| F66- 365651 | 200 | 40 | 5,0 | 4,0 | 1 | |
| F66- 365652 | 200 | 45 | 5,0 | 4,0 | 1 | |
| F66- 365653 | 200 | 50 | 5,0 | 4,0 | 1 | |
| F99- 366415 | 150 | 40 | 6,0 | 5,0 | 1 | |
| F99- 366416 | 150 | 45 | 6,0 | 5,0 | 1 | |
| F99- 366417 | 150 | 50 | 6,0 | 5,0 | 1 | |
| F99- 366645 | 175 | 40 | 6,0 | 5,0 | 1 | |
| F99- 366646 | 175 | 45 | 6,0 | 5,0 | 1 | |
| F99- 366647 | 175 | 50 | 6,0 | 5,0 | 1 | |
| F99- 366651 | 200 | 40 | 6,0 | 5,0 | 1 | |
| F99- 366652 | 200 | 45 | 6,0 | 5,0 | 1 | |
| F99- 366653 | 200 | 50 | 6,0 | 5,0 | 1 | |

SCHANZ SCREWS: SELF TAPPING, WITH HELICAL FLUTE




| Order Nr. | L= Length [mm] | A= Thread [mm] | Ø [mm] | S [mm] | AMOUNT | PRICE |
|-------------|-------------------|-------------------|-----------|-----------|--------|-------|
| F33- 8-1300 | 60 | 20 | 3,0 | ▲ | 1 | |
| F33- 8-1302 | 80 | 20 | 3,0 | ▲ | 1 | |
| F66- 8-1310 | 150 | 50 | 5,0 | ▲ | 1 | |
| F66- 8-1312 | 170 | 50 | 5,0 | ▲ | 1 | |
| F66- 8-1314 | 200 | 50 | 5,0 | ▲ | 1 | |
| F99- 8-1320 | 150 | 50 | 6,0 | ▲ | 1 | |
| F99- 8-1322 | 170 | 50 | 6,0 | ▲ | 1 | |
| F99- 8-1324 | 200 | 50 | 6,0 | ▲ | 1 | |

THREADED KIRSCHNER WIRES



| Order Nr. | L= Length [mm] | A= Thread [mm] | Ø [mm] | S [mm] | AMOUNT | PRICE |
|-------------|-------------------|-------------------|-----------|-----------|--------|-------|
| F22- 18086 | 160 | 15 | 1,2 | | 10 | |
| F22- 7-1310 | 150 | 5 | 1,6 | | 10 | |
| F22- 7-1312 | 150 | 15 | 1,6 | | 10 | |
| F22- 7-1314 | 150 | 15 | 2,0 | | 10 | |
| F22- 7-1316 | 200 | 15 | 2,0 | | 10 | |
| F22- 311201 | 150 | 5 | 1,2 | | 10 | |
| F22- 311601 | 150 | 15 | 1,6 | | 10 | |
| F22- 312001 | 150 | 15 | 2,0 | | 10 | |
| F33- 312501 | 150 | 15 | 2,5 | | 10 | |
| F33- 313001 | 150 | 15 | 3,0 | | 10 | |

KIRSCHER WIRES WITH TROCAR TIP



| Order Nr. | L= Length | A= Thread | Ø | S | AMOUNT | PRICE |
|-------------|-----------|-----------|-----|---|--------|-------|
| F66- 311531 | 310 | | 1,5 | | 10 | |
| F66- 311534 | 340 | | 1,5 | | 10 | |
| F99- 311537 | 370 | | 1,5 | | 10 | |
| F99- 311540 | 400 | | 1,5 | | 10 | |